

Stylo™

The new way of writing skincare



The Stylo™ tube rewrites your product's success story in a choice of three different applicators: the Stylo™ Liner – meticulous, the Stylo™ Pipette – precise, and the Stylo™ Bristle, which softly strokes the skin. Apply liquid and semi-solid skincare products easily and precisely. The aesthetic ring ensures correct discharge and protects the product. Consumers intuitively understand how to turn the ring thanks to the arrow pictogram. The cone-shaped cap protects the applicator. Start a new product story with Stylo™.

Applications for the Stylo™ Liner

- Anti-aging applications, e. g. wrinkle filler
- Make-up applications
- Creams for small blemishes and impurities

Applications for the Stylo™ Pipette

- Serum for pigmentation and aging spots
- Anti-aging formulations for eyes, lips and décolleté
- Serum for hair ends
- Cuticle serum/oil

- Brightening tooth gel
- Dietary supplement
- Anti-pruritic products

Applications for the Stylo™ Bristle

- Concealers/lotions
- Nail care
- Eyebrow care
- Acne cream
- Disinfectant solutions for piercings
- Scar care
- Mouth and teeth products

Advantages

- Precise application for different viscous products
- Easy-to-clean, hygienic applicators
- Intuitive ergonomic tubes with «pen» look
- Ring, cover and cap can be designed to customer specifications
- The ring ensures discharge and product protection

Specifications

MATERIAL		DIMENSIONS	
Tube	Polyfoil®, PE, COEX, Sugarcane, Picea	Diameter	19 mm
Liner	TPE	Content	10 – 20 ml
Pipette	PET-G		
Bristle	MBA™ microbristle applicator LDPE, pharma grade		
Ring	PP		
Cover	PP		
Cap	PET-G		

DECORATION

Offset print
Silk screen print
Hotfoil stamping
Varnish: glossy, satin,
matte

CERTIFICATIONS

Quality
Food grade
Bristle pharma grade